

78018-22 Masterflex® Ismatec® Reglo Digital Pump with MasterflexLive®, 4-Channel, 8-Roller; 115/230 VAC



Max. flow rate	228 ml/min
Height	157 mm (6.3")
Length	221 mm (8.7")
Width	147 mm (5.8")
Pump Head	Cartridge
Power	90 - 260 VAC; 50/60 Hz
Max. RPM	160 rpm
Housing material	ABS plastic and aluminium
No. of Channels	4
No. of rollers	8
Control Type	Digital
Max. flow rate per channel	57 ml/min
Min. Flow Rate	0.0002 ml/min
Drive IP rating	IP 31
Accuracy	±0.1%
Anti-drip	Yes
Cloud connected	Yes
Digital dispensing	Yes
Drive	Reglo
Duty Cycle	Continuous
Min. RPM	0.2 rpm
Open-head sensor	No
Remote I/O	RJ45 Ethernet; USB; 0 to 20 mA. 4 to 20 mA. 0 to 10 V. start/stop. direction. and prime via DB-25
Reversible motor	Yes
Series	Reglo
Tubing sizes accepted	Three-stop 0.13 mm ID to 3.17 mm ID
Washdown capable	No
Max. Pressure	1.5 bar (22 psi)
Min. flow rate per channel	0.0002 ml/min
Frequency	50/60 Hz
Voltage	90 - 260 VAC
Warranty	2 years

Distributed by:

Fiedler Scientific Instruments, s.r.o.
H. Malířové 1, 638 00 Brno, Czech Republic
Telephone: (+420) 541 222 795

Email: info@lab-eu.com info@pristroje.cz