

# Cole-Parmer®

## Disposable and Reusable Sampling Tools

- Sampling tools for free-flowing to viscous liquids
- Single-point, multi-level sampling for loose granules to cohesive powders
- Wide range of sampling tools commonly used in the lab

[coleparmer.com](http://coleparmer.com)



Distributed by:

Fiedler Scientific Instruments, s.r.o.  
H. Malířové 1, 638 00 Brno, Czech Republic  
Telephone: (+420) 541 222 795

Email: [info@lab-eu.com](mailto:info@lab-eu.com) [info@pristroje.cz](mailto:info@pristroje.cz)


## Disposable and Reusable Sampling Tools

### Benefits of Disposable Sampling Tools


- No risk of cross-contamination
- Available sterilized and individually packaged
- Most are manufactured in a cleanroom setting
- Lot-specific documentation available

### Benefits of Reusable Sampling Tools


- Includes long-lasting materials such as 316L stainless steel and PTFE
- Materials are easy to clean
- Certificates of Conformity available




## Liquid Sampling Tools Ordering Information

Description	Disposable/ Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Liquid Sample Collector, HDPE	Disposable	No	100 mL	500 mm	20/pk	06286-57	
	Disposable	Yes	100 mL	500 mm	20/pk	06286-58	
	Disposable	No	190 mL	1000 mm	20/pk	06286-59	
	Disposable	Yes	190 mL	1000 mm	20/pk	06286-60	


- Ideal for sampling liquids, thin oils, and ointments

Description	Disposable/ Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Viscous Sample Collector, HDPE	Disposable	No	190 mL	1000 mm	20/pk	06043-13	
	Disposable	Yes	190 mL	1000 mm	20/pk	06043-14	
	Disposable	No	100 mL	500 mm	20/pk	06043-11	
	Disposable	Yes	100 mL	500 mm	20/pk	06043-12	

- Large bore hole and oversized handle are ideal for viscous materials like honey, creams, pastes, and thick oils


Description	Disposable/ Reusable	Sterile	Capacity	Length x Diameter	Quantity	Item Number	
Ice Sampler, 316L Stainless Steel	Reusable	No	40 mL	125 mm x 20 mm	1/ea	06286-01	
	Reusable	No	400 mL	210 mm x 50 mm	1/ea	06286-02	

- Ideal for demanding applications such as sampling wax and frozen products
- Designed to fit most common commercial drills


Description	Disposable/ Reusable	Sterile	Capacity	Diameter	Quantity	Color	Item Number	
Ladle, PP	Disposable	Yes	50 mL	43 mm	40/pk	White	06286-34	
Ladle, HDPE	Disposable	Yes	50 mL	48 mm	40/pk	White	06286-35	
	Disposable	Yes	50 mL	48 mm	40/pk	Blue	06286-36	

- Easy to use with liquids, creams, and pastes


## Powder Sampling Tools Ordering Information

Description	Disposable/ Reusable	Sterile	Number of Openings	Capacity	Length	Quantity	Item Number	
Slot Sampler, 316L Stainless Steel	Reusable	No	3 openings	45 mL	300 mm	1/ea	06286-03	
	Reusable	No	3 openings	110 mL	600 mm	1/ea	06286-04	
	Reusable	No	4 openings	190 mL	1000 mm	1/ea	06286-05	
	Reusable	No	5 openings	300 mL	1500 mm	1/ea	06286-06	
	Reusable	No	5 openings	370 mL	2000 mm	1/ea	06286-07	
Slot Sampler, Aluminum	Reusable	No	3 openings	80 mL	600 mm	1/ea	06286-08	
	Reusable	No	5 openings	135 mL	1000 mm	1/ea	06286-09	
	Reusable	No	7 openings	190 mL	1500 mm	1/ea	06286-10	
	Reusable	No	7 openings	245 mL	2000 mm	1/ea	06286-11	
	Reusable	No	10 openings	350 mL	2500 mm	1/ea	06286-12	
Slot Sampler, HDPE	Disposable	Yes	1 opening	50 mL	620 mm	20/pk	06286-21	
	Disposable	Yes	1 opening	50 mL	970 mm	20/pk	06286-22	
	Disposable	Yes	3 openings	100 mL	620 mm	20/pk	06286-19	
	Disposable	Yes	3 openings	170 mL	970 mm	20/pk	06286-20	

- Ideal for free-flowing powder


Description	Disposable/ Reusable	Sterile	Number of Openings	Length	Quantity	Item Number	
Grain Spear, Aluminum	Reusable	No	4 openings	1500 mm	1/ea	06286-13	
	Reusable	No	8 openings	2000 mm	1/ea	06286-14	

- Ideal for food and environmental sampling
- Pockets positioned in spiral; bottom pocket opens and fills first, giving a more representative sample


Description	Disposable/ Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Sleeve Sampler Chamber, Stainless Steel	Reusable	No	100 mL	190 mm	1/ea	06286-15	
Sleeve Sampler Handle, Stainless Steel	Reusable	No	—	1000 mm	1/ea	06286-16	

- Use chamber and handle together for sampling at greater depths


Powder Sampling Tools Ordering Information (continued)


Description	Disposable/Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Powder Sample Collector, HDPE	Disposable	Yes	10 mL	500 mm	20/pk	06286-30	
	Disposable	Yes	10 mL	1000 mm	20/pk	06286-31	
	Disposable	Yes	20 mL	500 mm	20/pk	06286-32	
	Disposable	Yes	20 mL	1000 mm	20/pk	06286-33	

- Ideal for free-flowing powders

Description	Disposable/Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Powder Lance, HDPE	Disposable	No	200 mL	950 mm	20/pk	06286-41	
	Disposable	Yes	200 mL	950 mm	20/pk	06286-42	
	Disposable	No	75 mL	400 mm	20/pk	06286-43	
	Disposable	Yes	75 mL	400 mm	20/pk	06286-44	
	Disposable	No	125 mL	600 mm	20/pk	06286-45	
	Disposable	Yes	125 mL	600 mm	20/pk	06286-46	

- Ideal for cohesive powders




Description	Disposable/Reusable	Sterile	Capacity	Length	Quantity	Item Number	
Powder Spatula, PS	Disposable	Yes	150 mL	230 mm	100/pk	06264-19	
	Disposable	Yes	250 mL	330 mm	100/pk	06264-24	

Description	Disposable/Reusable	Dimensions	Item Number	
Cleaning Brush, Nylon Bristles	Reusable	1100 mm length, up to 30 mm diameter	06286-17	
	Reusable	1100 mm length, up to 45 mm diameter	06286-18	

- Ideal for cleaning stainless steel and plastic samplers

## General Lab Sampling Tools Ordering Information






### Scoops

Description	Disposable/ Reusable	Sterile	Capacity	Quantity	Item Number	
Industrial Scoop, Stainless Steel	Reusable	No	100 mL	1/ea	06700-63	
	Reusable	No	200 mL	1/ea	06700-64	
	Reusable	No	300 mL	1/ea	06700-65	
	Reusable	No	500 mL	1/ea	06700-66	
	Reusable	No	1000 mL	1/ea	06700-67	
	Reusable	No	1550 mL	1/ea	07210-00	
Pharma Scoop, 316L Stainless Steel	Reusable	No	10 mL	1/ea	07205-10	
	Reusable	No	50 mL	1/ea	07205-11	
	Reusable	No	100 mL	1/ea	07205-12	
	Reusable	No	250 mL	1/ea	07205-13	
	Reusable	No	500 mL	1/ea	07205-14	
	Reusable	No	750 mL	1/ea	07205-15	
	Reusable	No	1000 mL	1/ea	07205-16	
	Reusable	No	1250 mL	1/ea	07205-17	
	Reusable	No	1700 mL	1/ea	07205-18	
	Reusable	No	2500 mL	1/ea	07205-19	
	Reusable	No	5000 mL	1/ea	07205-20	
Round Scoop, Aluminum	Reusable	No	350 mL	1/ea	07207-02	
	Reusable	No	1700 mL	1/ea	07207-08	

- Pharma scoops feature hygienic design ideal for pharmaceutical industry; fully welded with no dead spaces or crevices and made from high quality corrosive-resistant stainless steel

General Lab Sampling Tools Ordering Information (continued)

Scoops

Description	Disposable/ Reusable	Sterile	Capacity	Quantity	Item Number	
Scoop, PP	Reusable	No	2 mL	12/pk	06287-27	
	Reusable	No	5 mL	12/pk	06287-28	
	Reusable	No	10 mL	12/pk	06287-29	
	Reusable	No	25 mL	12/pk	06287-30	
	Reusable	No	50 mL	12/pk	06287-31	
	Reusable	No	100 mL	12/pk	06287-32	
	Reusable	No	250 mL	12/pk	06287-33	
	Reusable	No	500 mL	12/pk	06287-34	
	Reusable	No	1000 mL	6/pk	06287-35	
Scoop, PTFE	Reusable	No	24 mL	1/ea	06300-07	
Long-Handled Scoop, HDPE	Disposable	No	50 mL	50/pk	06264-10	
	Disposable	No	100 mL	50/pk	06264-11	
Round-Bottom Scoop, HDPE	Reusable	No	250 mL	5/cs	06276-01	
	Reusable	No	500 mL	5/cs	06276-03	
	Reusable	No	1000 mL	5/cs	06276-04	
	Reusable	No	1500 mL	5/cs	06276-05	
Round Scoop, HDPE	Reusable	No	2000 mL	1/ea	06276-10	

- Polypropylene (PP) scoops are autoclavable
- Polystyrene (PS) sampling scoops are sterilized and designed with ergonomic handles

General Lab Sampling Tools Ordering Information (continued)

Scoops

Description	Disposable/ Reusable	Sterile	Capacity	Quantity	Item Number	
Scoop, PS	Disposable	Yes	30 mL	100/cs	06260-00	
	Disposable	Yes	60 mL	100/cs	06260-01	
	Disposable	Yes	120 mL	100/cs	06260-02	
	Disposable	Yes	250 mL	100/cs	06260-03	
	Disposable	Yes	500 mL	50/cs	06260-04	
	Disposable	Yes	1000 mL	25/cs	06260-05	
Scoop, PS, Double-bagged	Disposable	Yes	60 mL	50/cs	06302-00	
	Disposable	Yes	120 mL	50/cs	06302-01	
	Disposable	Yes	250 mL	25/cs	06302-02	
Scoop, PS, with lid	Disposable	Yes	60 mL	50/cs	06302-03	
	Disposable	Yes	120 mL	50/cs	06302-04	

- Polystyrene (PS) sampling scoops are gamma sterilized and FDA grade
- Ergonomic handles are curved at the bottom for better grip
- Flat bottoms provide stability on benchtop or scale
- Accurate volumes on all sizes



## General Lab Sampling Tools Ordering Information (continued)

### Spoons


Description	Disposable/ Reusable	Sterile	Capacity	Quantity	Item Number	
Volumetric Spoon, HDPE	Disposable	Yes	1 mL	100/pk	06264-00	
	Disposable	Yes	2.5 mL	100/pk	06264-01	
	Disposable	Yes	5 mL	100/pk	06264-02	
	Disposable	Yes	10 mL	100/pk	06264-03	
	Disposable	Yes	20 mL	100/pk	06264-04	
	Disposable	Yes	30 mL	100/pk	06264-05	
Spoon, PS	Disposable	Yes	2.5 mL	100/cs	06287-39	
	Disposable	Yes	10 mL	100/cs	06287-40	
Sampling spoon, PS	Disposable	Yes	1.23 mL (¼ tsp)	200/cs	06263-10	
	Disposable	Yes	2.46 mL (½ tsp)	200/cs	06263-20	
	Disposable	Yes	4.93 mL (1 tsp)	200/cs	06263-30	
	Disposable	Yes	14.79 mL (1 tbsp)	200/cs	06263-40	
Sampling spoon, PS, Double-bagged	Disposable	Yes	1.23 mL (¼ tsp)	50/cs	06302-13	
	Disposable	Yes	2.46 mL (½ tsp)	50/cs	06302-12	
	Disposable	Yes	4.93 mL (1 tsp)	50/cs	06302-11	
	Disposable	Yes	14.79 mL (1 tbsp)	50/cs	06302-10	
Long-Handled Spoon, HDPE	Disposable	No	10 mL	100/pk	06286-47	
	Disposable	Yes	10 mL	100/pk	06286-48	
Double-Ended Spoon, 304 Stainless Steel	Reusable	No	0.5 mL and 2 mL	1/ea	06700-70	
	Reusable	No	0.8 mL and 2.8 mL	1/ea	06700-71	
	Reusable	No	1.6 mL and 3 mL	1/ea	06700-72	
	Reusable	No	1.8 mL and 5 mL	1/ea	06700-73	
	Reusable	No	2 mL and 8 mL	1/ea	06700-74	

- Polystyrene (PS) sampling spoons are gamma sterilized and FDA grade
- Polystyrene (PS) sampling spoons feature ergonomic handles with curve for better grip
- Polystyrene (PS) sampling spoons have flat bottoms for stability on benchtop or scale

General Lab Sampling Tools Ordering Information (continued)



Spatulas

Description	Disposable/ Reusable	Sterile	Length	Quantity	Item Number	
Spatula Macro Spoon, PP	Reusable	No	150 mm	10/pk	06287-36	
	Reusable	No	180 mm	10/pk	06287-37	
	Reusable	No	210 mm	10/pk	06287-38	
Spatula Macro Spoon, Stainless Steel	Reusable	No	120 mm	3/pk	06287-04	
	Reusable	No	150 mm	3/pk	06287-05	
	Reusable	No	180 mm	3/pk	06287-06	
	Reusable	No	210 mm	3/pk	06287-07	
	Reusable	No	250 mm	3/pk	06287-08	
	Reusable	No	300 mm	3/pk	06287-09	
Spatula Micro Spoon, Stainless Steel	Reusable	No	150 mm	3/pk	06287-16	
	Reusable	No	180 mm	3/pk	06287-17	
	Reusable	No	210 mm	3/pk	06287-18	
Double-Ended Flat Spatula, Stainless Steel	Reusable	No	130 mm	3/pk	06287-10	
	Reusable	No	155 mm	3/pk	06287-11	
	Reusable	No	180 mm	3/pk	06287-12	
	Reusable	No	210 mm	3/pk	06287-13	
	Reusable	No	255 mm	3/pk	06287-14	
	Reusable	No	300 mm	3/pk	06287-15	

Description	Disposable/ Reusable	Sterile	Capacity	Blade Width	Quantity	Item Number	
Flat Spatula, PS	Disposable	Yes	—	13 mm	100/pk	06264-20	
Micro Spatula, PS	Disposable	Yes	2 mL	13 mm	100/pk	06264-21	
Narrow Micro Spatula, PS	Disposable	Yes	0.5 mL	6 mm	100/pk	06264-23	
Very-Narrow Micro Spatula, PS	Disposable	Yes	0.3 mL	4 mm	100/pk	06264-22	


## General Lab Sampling Tools Ordering Information (continued)

### Spatulas (continued)





Description	Disposable/ Reusable	Sterile	Length	Quantity	Item Number	
Lab Spatula, PTFE-Coated Stainless Steel	Reusable	No	13 mm spoon end, 38 mm flat round end	1/ea	06369-10	
Lab Spatula, PTFE-Coated Stainless Steel	Reusable	No	38 mm flat taper, 32 mm flat round end	1/ea	06369-11	

- PTFE coating prevents chemicals from sticking

### Scrapers

Description	Disposable/ Reusable	Sterile	Blade Width	Quantity	Item Number	
Scraper, HDPE	Disposable	No	75 mm	100/pk	06286-49	
	Disposable	Yes	75 mm	100/pk	06286-50	
	Disposable	No	110 mm	100/pk	06286-51	
	Disposable	Yes	110 mm	100/pk	06286-52	


### Forceps

Description	Disposable/ Reusable	Sterile	Length	Tip	Quantity	Item Number	
Forcep, Reinforced Nylon	Disposable	No	122 mm	Serrated, blunt	100/pk	06443-15	
	Disposable	No	111 mm	Serrated, pointed	100/pk	06443-25	
	Disposable	Yes	122 mm	Serrated, blunt	100/pk	06443-26	
Forcep, PS	Disposable	Yes	123 mm	Serrated, blunt	100/pk	06433-20	

- Non-sterile forceps are autoclavable; the sterile forceps are individually wrapped with easy to open peel-packs
- Choose serrated pointed tips for handling smaller objects.
- Forcep handles have ridges for secure grip.

General Lab Sampling Tools Ordering Information (continued)

Funnels

Description	Disposable/ Reusable	Sterile	Size	For use with	Stem OD	Quantity	Item Number	
Antistatic Weighing Funnel, PP	Disposable	No	Micro	96-well microtubes	6.21 mm	96/pk	06127-00	
	Disposable	No	Small	5, 10, 15, 20 mL flasks	7.4 mm	50/pk	06127-02	
	Disposable	No	Medium	100, 200, 250 mL flasks	10.35 mm	50/pk	06127-04	
	Disposable	No	Large	500, 1000 mL flasks	17 mm	20/pk	06127-06	
	Disposable	No	X-Large	1000, 2000 mL flasks	19.1 mm	10/pk	06127-08	

- Designed to fit volumetric flasks
- Autoclavable
- Feature flat bottoms for stability on benchtops, balances, or scales

coleparmer.com

USA:  
+1.800.323.4340  
+1.847.549.7600  
4890.1

Canada: +1.800.363.5900  
China: +86.21.5109.9909  
France: +33 (0) 1486 37800  
Germany: +49 (0) 9377 92030

India: +1.800.266.1244  
Italy: +39 (0) 1313 89513  
UK: +44 (0) 1480 272279  
All other countries: +1.847.549.7600

**Cole-Parmer®**  
an antylia scientific company

Find out more at [antylia.com](http://antylia.com)

Distributed by:

Fiedler Scientific Instruments, s.r.o.  
H. Malířové 1, 638 00 Brno, Czech Republic  
Telephone: (+420) 541 222 795

Email: [info@lab-eu.com](mailto:info@lab-eu.com) [info@pristroje.cz](mailto:info@pristroje.cz)